

K&H Coolin' Gel Pet Cot™

Care & Instructions

Made in China

K&H Manufacturing, LLC
Colorado Springs, CO
www.khmf.com | (719) 591-6950



Thank you for purchasing K&H Manufacturing's Coolin' Gel Pet Cot™

- 1** Fold out cover so it is flat.
- 2** Be sure all 4 frame tubes are threaded through the sleeves sewn along edges of the cover.
- 3** Insert all 4 frame tubes into the nylon corner supports.

With proper care, your K&H Coolin' Gel Pet Cot™ should provide years of comfort and service.

Care Instructions:
Wipe clean with a damp cloth.



K&H Limited One Year Warranty

Our warranty is as follows: K&H Manufacturing warrants to the original consumer, that this product will be free of defects and workmanship for a period of one year from the date of purchase. The manufacturer's liability hereunder is limited to replacement of the product. This warranty is void if the product has been damaged by accident, unreasonable use, neglect, tampering, or other causes not arising from defects in material or workmanship.

To exercise your warranty: If, during the warranty period, your product malfunctions, return it to the place of purchase with your **original purchase receipt**. If this cannot be done, return the product, postage prepaid, to K&H Manufacturing, with **dated proof of purchase or original date receipt**. A return authorization number, obtained by calling (719) 591-6950 must accompany all returns to manufacturer. Please include a note explaining the problem when returning the product.